

POTEE Polyclonal Antibody

Description

| | |
|--------------------------------|---|
| Product type | Primary Antibody |
| Code | BT-AP13237 |
| Host | Rabbit |
| Isotype | IgG |
| Size | 100ul, 50ul, 20ul |
| Immunogen | Synthesized peptide derived from human protein . at AA range: 350-430 |
| Mol wt | N/A |
| Species reactivity | Human, Rat, Mouse |
| Clonality | Polyclonal |
| Recommended application | WB, ELISA |
| Concentration | 1 mg/ml |
| Full name | POTE ankyrin domain family member E |
| Synonyms | POTE ankyrin domain family member E ;ANKRD26-like family C member 1A;Prostate, ovary, testis-expressed protein on chromosome 2;POTE-2 |

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Belongs to the actin family.,Contains 5 ANK repeats.,In the C-terminal section; belongs to the actin family.,In the N-terminal section; belongs to the POTE family.,

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com