

UN45A Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13376
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human UN45A
Mol wt	103840
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	UN45A
Synonyms	UN45A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


UN45A plays a role in cell proliferation and myoblast fusion| binds progesterone receptor (PGR; MIM 607311) and HSP90 (HSPCA; MIM 140571)| and acts as a regulator of the progesterone receptor chaperoning pathway (Price et al.| 2002

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com