

RPA1 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP13893
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	DNA-directed RNA polymerase I subunit RPA1
Synonyms	DNA-directed RNA polymerase I subunit RPA1 ;RNA polymerase I subunit A1;EC 2.7.7.6;A190;DNA-directed RNA polymerase I largest subunit;DNA-directed RNA polymerase I subunit A;RNA polymerase I

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The protein encoded by this gene is the largest subunit of the RNA polymerase I complex. The encoded protein represents the catalytic subunit of the complex, which transcribes DNA into ribosomal RNA precursors. Defects in this gene are a cause of the Cincinnati type of acrofacial dysostosis.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com