

C1T9B Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP14574
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from part region of human protein
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Complement C1q and tumor necrosis factor-related protein 9B
Synonyms	Complement C1q and tumor necrosis factor-related protein 9B ;C1q/TNF-related protein 9B;CTRP9B;Complement C1q and tumor necrosis factor-related protein 9-like

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com