

GPR97 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP15774
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human protein . at AA range: 210-290
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Probable G-protein coupled receptor 97
Synonyms	Probable G-protein coupled receptor 97 ;G-protein coupled receptor PGR26

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year