

Lyn Monoclonal Antibody

Description

Product type	Primary Antibody
Code	BT-MCA0857
Host	Mouse
Isotype	IgG
Size	50ul, 100ul
Immunogen	Purified recombinant fragment of Lyn expressed in E. Coli.
Mol wt	N/A
Species reactivity	Human
Clonality	Monoclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Tyrosine-protein kinase Lyn
Synonyms	LYN; JTK8; Tyrosine-protein kinase Lyn; Lck; Yes-related novel protein tyrosine kinase; V-yes-1 Yamaguchi sarcoma viral related oncogene homolog; p53Lyn; p56Lyn

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.


Recommended Dilution

ELISA: 1:10000

WB: 1:500 - 1:2000

Not yet tested in other applications.

Images



Western Blot analysis using Lyn Monoclonal antibody againsts HL60 (1) L540 (2) SLLP-M2 (3) SEM (4) and Ramos (5) cell lysate.

Storage

-20°C for one year