

KSHV ORF62 Monoclonal Antibody

Description

Product type	Antibody
Code	BT-MCA3374
Host	Mouse
Isotype	Mouse IgM
Size	100 μ L, 50 μ L
Immunogen	Purified recombinant fragment of human KSHV ORF62 expressed in E. Coli.
Mol wt	N/A
Species reactivity	KSHV
Clonality	Monoclonal
Recommended application	IHC
Concentration	N/A
Full name	N/A
Synonyms	ORF62

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Kaposi's sarcoma-associated herpesvirus (KSHV) belongs to the gamma-(2)-herpesvirus subfamily and has been closely linked to the Kaposi.

Recommended Dilution


WB: 1:500 - 1:2000

IHC-p: 1:200 - 1:1000

ELISA: 1:10000

Not yet tested in other applications.

Images



Immunocytochemistry analysis of TPA induced BCBL-1 cells (A) and uninduced BCBL-1 cells (B) using KSHV ORF62 mouse mAb with AEC staining.



Western blot analysis using KSHV ORF62 mouse mAb against KSHV ORF62 recombinant protein.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term.