



Optimize Your Research

Recombinant Cynomolgus B-cell Maturation Protein , C-Fc, Human Cells protein

Code CD00060

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

Source

Human Cells

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

Lyophilized from a 0.2 µm filtered solution of 50 mM Tris, 100 mM Glycine, pH 7.5.

Endotoxin

Less than 0.1 ng/µg (1 EU/µg) as determined by LAL test.

Reconstitution

Dissolve the lyophilized protein in distilled water.

If you have any question on order please contact us via: order@bt-laboratory.com; technical assistance please contact us via: support@bt-laboratory.com More product visit www.bt-laboratory.com

5/F 2 Bldg., 501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China
Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com