



Optimize Your Research

Recombinant Mouse Tissue Inhibitor of Metalloproteinases 2 (27-220), C-6His, Human Cells protein

Code CD00198

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

Source

Human Cells

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.

Endotoxin

Less than 0.1 ng/µg (1 EU/µg) as determined by LAL test.

If you have any question on order please contact us via: order@bt-laboratory.com; technical assistance please contact us via: support@bt-laboratory.com More product visit www.bt-laboratory.com

5/F 2 Bldg., 501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China
Tel: 86 21 31007137 | E-mail: support@bt-laboratory.com | www.bt-laboratory.com