



Optimize Your Research

Recombinant Mouse Ciliary Neurotrophic Factor protein

Code CD00514

Storage: This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

Intended Use

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

Source

Escherichia coli.

Molecular Weight

Approximately 22.5 kDa, a single non-glycosylated polypeptide chain containing 197 amino acids.

Purity

> 95 % by SDS-PAGE and HPLC analyses.

Biological Activity

Fully biologically active when compared to standard. The ED50 as determined by a cell proliferation assay using human TF-1 cells is less than 30 ng/ml, corresponding to a specific activity of $> 3.3 \times 10^4$ IU/mg.

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

Lyophilized from a 0.2 µm filtered concentrated solution in 2 × PBS, pH 7.4, 2 % trehalose.

Sequence

AFAEQSPLTL HRRDLCSRSI WLARKIRSDL TALMESYVKH QGLNKNISLD SVDGVPVAST DRWSEMTEAE RLQENLQAYR TFQGMLTKLL
EDQRVHFTPT EGDHFHQAIHT LTLQVSAFAY QLEELMALLE QKVPEKEADG MPVTIGDGGGL FEKKLWGLKV LQELSQWTVR SIHDLRVISS
HHMGISAHES HYGAKQM

Endotoxin

Less than 1 EU/µg of rMuCNTF as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and stored at ≤ -20 °C. Further dilutions should be made in appropriate buffered solutions.

If you have any question on order please contact us via: order@bt-laboratory.com; technical assistance please contact us via: support@bt-laboratory.com More product visit www.bt-laboratory.com