



*Optimize Your Research*

# **Recombinant Human Epidermal Growth Factor Receptor/Receptor Tyrosine Protein Kinase ErbB1 , C-6His, Human Cells protein**

**Code CD01659**

**Storage:** This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2-8°C. For maximal stability, apportion the reconstituted preparation into working aliquots and store at -20°C to -70°C. Avoid repeated freeze/thaw cycles.

## **Intended Use**

This product is for research use only, not for use in diagnosis procedures. It is highly recommended to read this instruction entirely before the use.

## **Source**

Human Cells

## **Purity**

Greater than 95% as determined by reducing SDS-PAGE.

## **Physical Appearance**

Sterile Filtered White lyophilized (freeze-dried) powder.

## **Formulation**

Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.

## **Endotoxin**

Less than 0.1 ng/µg (1 EU/µg) as determined by LAL test.

## **Reconstitution**

Dissolve the lyophilized protein in distilled water.

If you have any question on order please contact us via: [order@bt-laboratory.com](mailto:order@bt-laboratory.com); technical assistance please contact us via: [support@bt-laboratory.com](mailto:support@bt-laboratory.com) More product visit [www.bt-laboratory.com](http://www.bt-laboratory.com)

5/F 2 Bldg., 501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China  
Tel: 86 21 31007137 | E-mail: [support@bt-laboratory.com](mailto:support@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)