



# Try Our One-Step ELISA Kit

Based on sandwich ELISA assay BT LAB's developed One-Step ELISA assay to be more simple. Unlike conventional sandwich ELISA, One-Step ELISA kits generate data within just 90 minutes. Without diluting, add the sample directly into the well, just in one incubation and washing step, the complete sandwich complex forms in the well and binds to the plate.

## Feaures of One-Step Kit

- Generate data in 90 mins
- Ready to use with No sample dilution
- Antibody pair for non-specific binding
- Validated in various biological samples

## Recommend Analytes

Rat Interleukin 10  
Rat Interleukin 6  
Rat TNF-alpha  
Mouse Interleukin 10  
Mouse Interleukin 6  
Human TNF-alpha

## Order Now !

86 (021) 3100 7137

[order@bt-laboratory.com](mailto:order@bt-laboratory.com)

visit us: [www.bt-laboratory.com](http://www.bt-laboratory.com)