

QUANTITATIVE ELISAs

a wide selection for measuring molecules in biological fluids

Cytokine

Adiponectin, BDNF, BLC, Beta-NGF, Eotaxin-2, E-Selectin, EGFR, FGF2, GDNF, GM-CSF, GRO, IL-1 alpha, IL-1 beta, IL-2, IL-4, IL-5, IL-6, IL-7, IL-8, IL-10, IL-13, IL-17, IFN-gamma, ICAM-1, IGFBP-1, IGFBP-2, IGFBP-3, HGF, L-Selectin, M-CSF, MCP-1, MMP-9, MIP-3 beta, RANTES, TGF-1, TNF-alpha, TNF-beta, VEGF-A

Chemokine

Axl, CTACK/CCL27, CCL1, CXCL16, ENA-78, Eotaxin-3, GCP-2, GRO, HCC-1, IL-8, IL-9, IL-17F, IL-28A, IL-29, IL-31, IP-10, I-TAC, LIF, Lymphotactin, MCP-2, MCP-3, MCP-4, MDC, MIF, MIP-3 alpha, MIP-3 beta, MPIF-2, NAP-2, Osteopontin, PARC, PF4, TARC, TECK, TSLP

Inflammation

ANP, AVP, BLC, CRP, Eotaxin-1, Eotaxin-2, ET-1, G-CSF, GM-CSF, ICAM-1, IFN-gamma, IL-1 alpha, IL-1 beta, IL-2, IL-4, IL-5, IL-6, IL-7, IL-8, IL-10, IL-11, IL-13, IL-15, IL-16, IL-17, MCP-1, M-CSF, MIG, MIP-1 alpha, MIP-1 beta, MPO, OSM, PDGF-BB, RANTES, SAA, TIMP-1, TIMP-2, TNF-alpha

Apoptosis

AKT, Annexin-5, APC, BAX, BCL-2, Caspase, Catalase, CD40, CD95, CD27, IGF1, IGF1R, GAPDH, HDAC1, HSP70, HIF1-alpha, PARP-1, PD-1

Growth Factor

Amphiregulin, ANGPT-1, BDNF, bFGF, BMP-2, BMP-4, BMP-5, BMP-7, Beta-NGF, CSF-3, EGF, EGFR, EG-VEGF, EPO, FGF-2, FGF-4, FGF-7, Fibronectin, GDF-15, GDNF, Growth Hormone, HB-EGF, HGF, IGF-1, IGF-2, IGFBP-1, IGFBP-2, IGFBP-3, IGFBP-4, IGFBP-6, IGF-1, Insulin, NGF R, NT-3, NT-4, Osteoprotegerin, PDGF-AA, PIGF, SCF, TGF-alpha, TGF-beta1, TGF-beta3, VEGFR2, VEGF-A, VEGF-D

Hormone

Adrenaline, Amylin, ACTH, Cortisol, Dopamine, Estradiol, Follistatin, FSH, LH, Progesterone, Prolactin, Resistin, Testosterone, TSH

Tumor Marker

EPCAM, HER2, HSP70, MTA1, RBP, RBP4, OSM, PAI1

Cardiac Marker

LDL, LDLR, Lipoprotein, LOX-1, Myoglobin, MYL1, Troponin T, TNNI2, TNNI3, Troponin I

Human Obesity

Adiponectin, AgRP, ANGPTL4, BDNF, Chemerin, CRP, GH, IFN-gamma, IGFBP-1, IGFBP-2, IGF-I, IL-10, IL-1 beta, IL-1 RA, IL-6, IL-8, Insulin, Leptin, VEGF, Lipocalin-2, OPG, PAI-1, PDGF-BB, Pepsinogen 1, Pepsinogen 2, Procalcitonin, Prolactin, RANTES, RBP4, Resistin, SAA, TGF-beta 1, TSP-1, TNF-alpha,

Metabolic

Adiponectin, Amylin, C-Peptide, FGF-21, FGF-23, Ghrelin, GLP-1, Glucagon, Proinsulin, PYY, Resistin, SAA-3

More ELISA panels, please visit us at www.bt-laboratory.com

BioassayTech Lab OPTIMIZE YOUR RESEARCH

Shanghai Korain Biotech Co. provides a wide range of enzyme-linked immunosorbent assay (ELISA) kits for the detection of thousands of different proteins and molecules including cytokines, growth factors, markers for infectious diseases, diabetes and tumor, drugs and small molecules.

- Broad dynamic detection range
- High specificity with no cross-reactivity
- Rigorous validation testing
- Technical resources available

Can't find the kit you need? Contact us at save@bt-laboratory.com to customize an affordable one for you! Choose from our catalog of 15 species to create an ELISA kit specific to your research.