

## Dylight 680, Goat Anti-Rabbit IgG

### Description

<b>Product type</b>	Secondary Antibody
<b>Code</b>	AS0414
<b>Host</b>	Goat
<b>Isotype</b>	N/A
<b>Size</b>	100ul, 1ml
<b>Immunogen</b>	N/A
<b>Mol wt</b>	N/A
<b>Species reactivity</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	IF, ICC, FCM, WB
<b>Concentration</b>	1mg/mL
<b>Full name</b>	Dylight 680, Goat Anti-Rabbit IgG
<b>Synonyms</b>	N/A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

Secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

### Recommended Dilution

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)