

Goat Anti-Rabbit IgG (H&L) FITC

Description

Product type Secondary Antibody

Code AS0472

Host Goat

Isotype N/A

Size 100ul, 1ml

Immunogen N/A

Mol wt N/A

Species reactivity Rabbit

Clonality Polyclonal

Recommended application IF, FCM

Concentration 1 mg/mL

Full name Goat Anti-Rabbit IgG (H&L) FITC

Synonyms

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

Recommended Dilution

IF: 1: 250 - 1: 2000 FCM: 1: 250 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year