

## Goat Anti-Mouse IgG (H&L)

### Description

<b>Product type</b>	Secondary Antibody
<b>Code</b>	AS0784
<b>Host</b>	Goat
<b>Isotype</b>	N/A
<b>Size</b>	100ul, 1ml
<b>Immunogen</b>	Other Immunogen Type corresponding to Mouse IgG
<b>Mol wt</b>	N/A
<b>Species reactivity</b>	Mouse
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	ELISA, IF, ICC, IHC
<b>Concentration</b>	4 mg/ml
<b>Full name</b>	Goat Anti-Mouse IgG (H&L)

### Synonyms

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

Immunoway secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

### Recommended Dilution

IHC: 1: 500 - 1: 1000

IF: 1: 5000 - 1: 10000

WB: 1: 5000 - 1: 20000

ELISA: 1: 10000 - 1: 20000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for one year