

Goat Anti-Mouse IgG (H&L) FITC

Description

Product type	Secondary Antibody
Code	AS0816
Host	Goat
Isotype	IgG1
Size	100ul, 1ml
Immunogen	N/A
Mol wt	N/A
Species reactivity	Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA, IF, IHC
Concentration	2mg/mL
Full name	Goat Anti-Mouse IgG (H&L) FITC
Synonyms	N/A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

IHVSecondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

Recommended Dilution

IF: 1: 5000 - 1: 10000

IHC: 1: 5000 - 1: 10000

WB: 1: 5000 - 1: 20000

ELISA: 1: 10000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year