

ALKBH5 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00070
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human ALKBH5
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	IHC-p, IF, WB
Concentration	1 mg/ml
Full name	ALKBH5
Synonyms	ALKBH5; ABH5; OFOX D; OFOX D1; Probable alpha-ketoglutarate-dependent dioxygenase ABH5; EC 1.14.11.-; Alkylated DNA repair protein alkB homolog 5

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Belongs to the alkB family.]

Recommended Dilution

WB: 1: 500 - 1: 2000

IHC-p: 1: 50 - 1: 200

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com