

ABHD14A Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00118
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD14A. AA range:201-250
Mol wt	29765
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ABHD14A Antibody
Synonyms	ABHD14A; Abhydrolase domain-containing protein 14A

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Possible role in granule neuron development.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 40000

Not yet tested in other applications.

Images



Western blot analysis of lysates from RAW264.7 cells, using ABHD14A Antibody. The lane on the right is blocked with the synthesized peptide.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com