

Actin beta Rabbit Polyclonal Antibody(Drosophilidae)

Description

Product type	Primary Antibody
Code	BT-AP00234
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Purified recombinant Full length Actin beta (Drosophilidae)
Mol wt	N/A
Species reactivity	Drosophilidae
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Actin beta
Synonyms	Actin beta; Beta Actin; β Actin

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


Recommended Dilution

WB: 1: 5000 - 1: 20000

ELISA: 1: 20000

Not yet tested in other applications.

Images



Western Blot analysis of Drosophilidae cells using Actin beta rabbit pAb(Drosophilidae). Secondary antibody was diluted at 1:20000

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com