

ADAM 12 (Cleaved-Glu208) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00258
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human ADAM 12 (Cleaved-Glu208)
Mol wt	99990
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ADAM 12
Synonyms	ADAM 12; Disintegrin and metalloproteinase domain-containing protein 12; ADAM 12; EC 3.4.24.-; Meltrin-alpha

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com