

ALG3 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00357
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human ALG3
Mol wt	48180
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	ALG3
Synonyms	Dol-P-Man:Man;5GlcNAc;2-PP-Dol alpha-1; 3-mannosyltransferase; EC 2.4.1.258;Asparagine-linked glycosylation protein 3 homolog;Dol-P-Man-dependent alpha;1-3-mannosyltransferase;Dolichyl-P-Man:M

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiading, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com