

Amphiregulin Polyclonal Antibody

Description

Product type Primary Antibody

Code BT-AP00413

Host Rabbit

Isotype IgG

Size 20ul, 50ul, 100ul

Immunogen The antiserum was produced against synthesized peptide derived from the Internal region of human

AREG/AREGB. AA range:91-140

Mol wt 27895

Species reactivity Human

Clonality Polyclonal

Recommended application WB, ELISA

Concentration 1 mg/ml

Full name Amphiregulin Antibody

Synonyms AREG; SDGF; AREGB; Amphiregulin; AR; Colorectum cell-derived growth factor; CRDGF

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


The amphiregulin encoded by AREG is a member of the epidermal growth factor family. It is an autocrine growth factor as well as a mitogen for astrocytes, Schwann cells and fibroblasts. It is related to epidermal growth factor (EGF) and transforming growth factor alpha (TGF-alpha). The protein interacts with the EGF/TGF-alpha receptor to promote the growth of normal epithelial cells, and it inhibits the growth of certain aggressive carcinoma cell lines. It also functions in mammary gland, oocyte and bone tissue development. This gene is associated with a psoriasis-like skin phenotype, and is also associated with other pathological disorders, including various types of cancers and inflammatory conditions.

Recommended Dilution

WB: 1: 500 - 1: 2000 ELISA: 1: 10000

Not yet tested in other applications.

Images



Western Blot analysis of HBE cells using Amphiregulin Polyclonal Antibody. Secondary antibody was diluted at

1:20000

Storage

-20°C for one year