

ANKRD20 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00451
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human ANKRD20A3. AA range:551-600
Mol wt	94108
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	ANKRD20 Antibody
Synonyms	ANKRD20A3; Ankyrin repeat domain-containing protein 20A3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

ANKRD20A2 (Ankyrin Repeat Domain 20 Family Member A2) is a Protein Coding gene. An important paralog of this gene is ANKRD18A.


Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 10000

Not yet tested in other applications.

Images



Immunohistochemistry analysis of paraffin-embedded human brain, using ANKRD20A3 Antibody. The picture on the right is blocked with the synthesized peptide.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com