

ApoF Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00536
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human APOF. AA range:221-270
Mol wt	33463
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ApoF Antibody
Synonyms	APOF; Apolipoprotein F; Apo-F; Lipid transfer inhibitor protein; LTIP

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

The product of APOF is one of the minor apolipoproteins found in plasma. This protein (apolipoprotein F) forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol.


Recommended Dilution

WB: 1: 500 - 1: 2000


ELISA: 1: 40000

Not yet tested in other applications.


Images




Western Blot analysis of A549 cells using ApoF Polyclonal Antibody diluted at 1:1000



Western blot analysis of lysates from HeLa, Jurkat, COLO, and A549 cells, using APOF Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using ApoF Polyclonal Antibody diluted at 1:1000



Western blot analysis of the lysates from COLO205 cells using APOF antibody.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com