

ApoL6 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00543
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from ApoL6 . at AA range: 240-320
Mol wt	38128
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ApoL6 Antibody
Synonyms	APOL6; Apolipoprotein L6; Apolipoprotein L-VI; ApoL-VI

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

APOL6 is a member of the apolipoprotein L gene family. The apolipoprotein L6 is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles.


Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 10000

Not yet tested in other applications.

Images



Storage

-20°C for one year