

ARHGEF16 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00604
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from ARHGEF16 . at AA range: 140-220
Mol wt	80105
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ARHGEF16 Antibody
Synonyms	ARHGEF16; EPHEXIN4; NBR; Rho guanine nucleotide exchange factor 16; Ephexin-4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Although the specific function of Rho guanine nucleotide exchange factor 16 is not known yet, it is thought to be involved in protein-protein and protein-lipid interactions.

Recommended Dilution

WB: 1: 500 - 1: 2000

ELISA: 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com