

BACH1/BRIP1(Phospho-Thr1133) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00715
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human BACH1 and BRIP1 (Thr1133)
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	BACH1/BRIP1
Synonyms	BACH1/BRIP1 ;Thr1133; Fanconi anemia group J protein; Protein FACJ; EC 3.6.4.13; ATP-dependent RNA helicase BRIP1; BRCA1-associated C-terminal helicase 1; BRCA1-interacting protein C-terminal helicase 1; BRCA1-interacting protein 1

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

N/A

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com