

Acrogranin Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00849
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Acrogranin
Mol wt	65230
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Acrogranin
Synonyms	Acrogranin; Granulins; Proepithelin; PEPI; Acrogranin; Paragranulin; Granulin-1; Granulin G; Granulin-2; Granulin F; Granulin-3; Granulin B; Granulin-4; Granulin A; Granulin-5; Granulin C; Granulin-6; Granulin D; Granulin-7; Granulin E;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:Defects in GRN are the cause of ubiquitin-positive frontotemporal dementia (UP-FTD)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com