

BCAS4 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00854
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human BCAS4. AA range:31-80
Mol wt	22758
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, IHC-p, ELISA
Concentration	1 mg/ml
Full name	BCAS4 Antibody
Synonyms	BCAS4; Breast carcinoma-amplified sequence 4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

BCAS4 (Breast Carcinoma Amplified Sequence 4) is a Protein Coding gene. An important paralog of this gene is BLOC1S4.

Recommended Dilution


WB: 1: 500 - 1: 2000

IHC: 1: 100 - 1: 300


ELISA: 1: 20000

Not yet tested in other applications.

Images



Immunohistochemistry analysis of paraffin-embedded human tonsil tissue, using BCAS4 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from K562 cells, using BCAS4 Antibody. The lane on the right is blocked with the synthesized peptide.

Storage

-20°C for one year