

C-Jun (Hydroxylated-p244) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP00894
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human C-Jun (Hydroxylated-p244)
Mol wt	36410
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	C-Jun
Synonyms	C-Jun ;Hydroxylated-p244; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Transcription factor that recognizes and binds to the enhancer heptamer motif 5'-TGA

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com