

## C16orf44 Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP01026
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	20ul, 50ul, 100ul
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human C16orf44. AA range:567-616
<b>Mol wt</b>	69896
<b>Species reactivity</b>	Human
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	IHC-p, ELISA
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	C16orf44 Antibody
<b>Synonyms</b>	KLHL36; C16orf44; Kelch-like protein 36

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

KLHL36 (Kelch Like Family Member 36) is a Protein Coding gene. An important paralog of this gene is KLHL22.

### Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for one year