

CaBP7 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP01079
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human CABP7. AA range:1-50
Mol wt	24453
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, IHC-p, ELISA
Concentration	1 mg/ml
Full name	CaBP7 Antibody
Synonyms	CABP7; CALN2; Calcium-binding protein 7; CaBP7; Calneuron II; Calneuron-2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Negatively regulates Golgi-to-plasma membrane trafficking by interacting with PI4KB and inhibiting its activity.

Recommended Dilution


WB: 1: 500 - 1: 2000

IHC: 1: 100 - 1: 300

ELISA: 1: 10000

Not yet tested in other applications.

Images



Western Blot analysis of COLO cells using CaBP7 Polyclonal Antibody diluted at 1:1000

Storage

-20°C for one year