

CDC20(Phospho-Ser51) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP01315
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human CDC20 (Ser51)
Mol wt	N/A
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	CDC20
Synonyms	CDC20 ;Ser51; Cell division cycle protein 20 homolog; p55CDC

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes| nuclear movement prior to anaphase and chromosome separation.

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com