

## CHM4B Rabbit Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP01477
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	100ul, 50ul, 20ul
<b>Immunogen</b>	Synthesized peptide derived from human CHM4B
<b>Mol wt</b>	24640
<b>Species reactivity</b>	Human, Mouse
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	Charged multivesicular body protein 4b
<b>Synonyms</b>	Charged multivesicular body protein 4b ;Chromatin-modifying protein 4b;CHMP4b;SNF7 homolog associated with Alix 1;SNF7-2;hSnf7-2;Vacuolar protein sorting-associated protein 32-2;Vps32-2;

**This product is for research use only, not for use in human, therapeutic or diagnostic procedure.**

### Background

#### Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

### Images

No images.

### Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhui Dist, Jiading, Zhejiang, China

Tel: 86 21 31007137 | E-mail: [save@bt-laboratory.com](mailto:save@bt-laboratory.com) | [www.bt-laboratory.com](http://www.bt-laboratory.com)