

Cdc50C Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP01599
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human TMEM30C. AA range:51-100
Mol wt	13039
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IF, ELISA
Concentration	1 mg/ml
Full name	Cdc50C Antibody
Synonyms	TMEM30C; CDC50C; Cell cycle control protein 50C; Transmembrane protein 30C

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

TMEM30C, also named as CDC50C, belongs to the CDC50/LEM3 family. TMEM30C is specifically expressed in testis.


Recommended Dilution

IF: 1: 200 - 1: 1000

ELISA: 1: 20000

Not yet tested in other applications.

Images



Immunofluorescence analysis of MCF7 cells, using TMEM30C Antibody. The picture on the right is blocked with the synthesized peptide.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com