

ALK(Phospho-Tyr1078) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP01622
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized phospho peptide around human ALK (Tyr1078)
Mol wt	N/A
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	Alkbh8
Synonyms	ALK ;Tyr1078 ; ALK tyrosine kinase receptor; EC 2.7.10.1; Anaplastic lymphoma kinase; CD antigen CD246

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes a receptor tyrosine kinase which belongs to the insulin receptor superfamily. This protein comprises an extracellular domain| an hydrophobic stretch corresponding to a single pass transmembrane region| and an intracellular kinase domain. It plays an important role in the development of the brain and exerts its effects on specific neurons in the nervous system. This gene has been found to be rearranged| mutated| or amplified in a series of tumours including anaplastic large cell lymphomas| neuroblastoma| and non-small cell lung cancer. The chromosomal rearrangements are the most common genetic alterations in this gene| which result in creation of multiple fusion genes in tumorigenesis| including ALK (chromosome 2)|EML4 (chromosome 2)| ALK/RANBP2 (chromosome 2)| ALK/ATIC (chromosome 2)| ALK/TFG (chromosome 3)| ALK/NPM1 (chromosome 5)| ALK/SQSTM1 (chromosome

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com