

ANGPTL3 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP01749
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human ANGPTL3
Mol wt	50600
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	ANGPTL3
Synonyms	ANGPTL3; Angiotensin-converting enzyme 2-related protein 3; Angiotensin-converting enzyme 2-related protein 3

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Contains 1 fibrinogen C-terminal domain. |tissue specificity: Expressed principally in liver. Weakly expressed in kidney. |

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com