

Collagen II alpha 1 (Cleaved-Ala1241) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP01770
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Collagen II α 1 (Cleaved-Ala1241)
Mol wt	163570
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Collagen II alpha 1
Synonyms	Collagen II α 1 ;Cleaved-Ala1241; Collagen alpha-1;II chain; Alpha-1 type II collagen; Collagen alpha-1;II chain; Chondrocalcin;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:Defects in COL2A1 are a cause of primary avascular necrosis of femoral head (ANFH)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com