

AP4M1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02180
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human AP4M1
Mol wt	49830
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	AP-4 complex subunit mu-1
Synonyms	AP4M1; MU ARP2; Ap4m4; MUARP2; SPG50; AP-4 complex subunit mu-1 ;AP-4 adapter complex mu subunit;Adapter-related protein complex 4 mu-1 subunit;Mu subunit of AP-4;Mu-adaptin-related protein 2;mu-ARP2;Mu4-adaptin;mu4

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Recommended Dilution

WB: 1: 1000 - 1: 2000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com