

CPXM2 Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02196
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	The antiserum was produced against synthesized peptide derived from human CPXM2. AA range:450-500
Mol wt	85898
Species reactivity	Human
Clonality	Polyclonal
Recommended application	IHC-p, ELISA
Concentration	1 mg/ml
Full name	CPXM2 Antibody
Synonyms	CPXM2; Inactive carboxypeptidase-like protein X2

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

May be involved in cell-cell interactions.

Recommended Dilution

IHC: 1: 100 - 1: 300

ELISA: 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com