

DC-LAMP Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02522
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthetic peptide from human protein at AA range: 261-310
Mol wt	26183
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, IHC-p, ELISA
Concentration	1 mg/ml
Full name	DC-LAMP Antibody
Synonyms	Lysosome-associated membrane glycoprotein 3 (LAMP-3) (Lysosomal-associated membrane protein 3) (DC-lysosome-associated membrane glycoprotein) (DC LAMP) (Protein TSC403) (CD antigen CD208)

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

CD302 is a C-type lectin receptor involved in cell adhesion and migration, as well as endocytosis and phagocytosis (Kato et al. 2007 [PubMed: 17947679]).

Recommended Dilution

WB: 1: 500 - 2000

IHC-p: 1: 500 - 200

ELISA: 1: 10000 - 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for one year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com