

F90A9 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02557
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from human F90A9
Mol wt	51040
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	F90A9
Synonyms	F90A9

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


FAM90A9 belongs to subfamily II of the primate-specific FAM90A gene family| which originated from multiple duplications and rearrangements (Bosch et al. 2007)

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4°C overnight

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com