

Fabp1 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02618
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human Fabp1
Mol wt	13970
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	Fabp1
Synonyms	Fabp1; Fatty acid-binding protein; liver; Fatty acid-binding protein 1; Liver-type fatty acid-binding protein; L-FABP

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

This gene encodes the fatty acid binding protein found in liver. Fatty acid binding proteins are a family of small| highly conserved| cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. This protein and FABP6 (the ileal fatty acid binding protein) are also able to bind bile acids. It is thought that FABPs roles include fatty acid uptake| transport| and metabolism.

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com