

FoxO1A (Acetyl Lys248) Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02719
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human FoxO1A (Acetyl Lys248)
Mol wt	72050
Species reactivity	Human, Mouse, Rat
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	FoxO1A
Synonyms	FoxO1A ;Acetyl Lys248; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

disease:Chromosomal aberrations involving FOXO1 are a cause of rhabdomyosarcoma 2 (RMS2)

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com