

GCP-2 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02801
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human GCP-2
Mol wt	12540
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	GCP-2
Synonyms	GCP-2; C-X-C motif chemokine 6; Chemokine alpha 3; CKA-3; Granulocyte chemotactic protein 2; GCP-2; Small-inducible cytokine B6; Small-inducible cytokine B6; N-processed variant 1; Small-inducible cytokine B6; N-processed variant 2; Small-inducible cytokine B6; N-processed variant 3;

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Chemotactic for neutrophil granulocytes. |online information: CXCL6 entry| Belongs to the intercrine alpha (chemokine Cx) family. |

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com