

GMCLL Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP02852
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human GMCLL
Mol wt	57860
Species reactivity	Human, Mouse
Clonality	Polyclonal
Recommended application	WB
Concentration	1 mg/ml
Full name	GMCLL
Synonyms	GMCLL

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background


This locus shares a high degree of identity with the multi-exon germ cell-less gene on chromosome 2. Despite its single-exon nature| this chromosome 5 locus contains an open reading frame that could putatively encode a full-length germ cell-less related protein.

Recommended Dilution

WB: 1: 500 - 1: 2000

Not yet tested in other applications.

Images



Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com