

IL-9 Rabbit Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03085
Host	Rabbit
Isotype	IgG
Size	100ul, 50ul, 20ul
Immunogen	Synthesized peptide derived from human IL-9
Mol wt	15840
Species reactivity	Human, Rat, Mouse
Clonality	Polyclonal
Recommended application	WB, ELISA
Concentration	1 mg/ml
Full name	IL-9
Synonyms	IL-9; Interleukin-9; IL-9; Cytokine P40; T-cell growth factor P40

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

Supports IL-2 independent and IL-4 independent growth of helper T-cells.[online information:Interleukin-9 entry|Belongs to the IL-7/IL-9 family.]

Recommended Dilution

WB: 1: 1000 - 1: 2000

ELISA: 1: 5000 - 1: 20000

Not yet tested in other applications.

Images

No images.

Storage

-20°C for 1 year

501 Changsheng S Rd, Nanhu Dist, Jiaxing, Zhejiang, China

Tel: 86 21 31007137 | E-mail: save@bt-laboratory.com | www.bt-laboratory.com