

FAM3D Polyclonal Antibody

Description

Product type	Primary Antibody
Code	BT-AP03151
Host	Rabbit
Isotype	IgG
Size	20ul, 50ul, 100ul
Immunogen	Synthesized peptide derived from the C-terminal region of human FAM3D.
Mol wt	24963
Species reactivity	Human
Clonality	Polyclonal
Recommended application	WB, IHC-p, ELISA
Concentration	1 mg/ml
Full name	FAM3D Antibody
Synonyms	FAM3D; Protein FAM3D

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

Background

FAM3D (Family With Sequence Similarity 3 Member D) is a Protein Coding gene. Diseases associated with FAM3D include narcolepsy. GO annotations related to this gene include cytokine activity. An important paralog of this gene is FAM3B.

Recommended Dilution


WB: 1: 500 - 1: 2000

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°over night

Storage

-20°C for one year