

## FAM3D Polyclonal Antibody

### Description

<b>Product type</b>	Primary Antibody
<b>Code</b>	BT-AP03151
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Size</b>	20ul, 50ul, 100ul
<b>Immunogen</b>	Synthesized peptide derived from the C-terminal region of human FAM3D.
<b>Mol wt</b>	24963
<b>Species reactivity</b>	Human
<b>Clonality</b>	Polyclonal
<b>Recommended application</b>	WB, IHC-p, ELISA
<b>Concentration</b>	1 mg/ml
<b>Full name</b>	FAM3D Antibody
<b>Synonyms</b>	FAM3D; Protein FAM3D

This product is for research use only, not for use in human, therapeutic or diagnostic procedure.

### Background

FAM3D (Family With Sequence Similarity 3 Member D) is a Protein Coding gene. Diseases associated with FAM3D include narcolepsy. GO annotations related to this gene include cytokine activity. An important paralog of this gene is FAM3B.

### Recommended Dilution


WB: 1: 500 - 1: 2000

IHC: 1: 100 - 1: 300

ELISA: 1: 40000

Not yet tested in other applications.

### Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night

### Storage

-20°C for one year